



# PLASTIC COMPONENTS FOR MEDICAL DEVICES

[www.bonini.com](http://www.bonini.com)







## PLASTIC COMPONENTS FOR MEDICAL DEVICES

### **Bag connectors**

A wide range of bag connectors for all market needs.



### **Assembled chambers**

A wide variety of assembled chambers with 1-2 way small and the solution filter



### **Silicone Products**

Our silicone medical products. Bonini quality serves innovation




### **Components**

All our components designed to meet your needs



# Assembled chambers



## GM21Fgap

Assembled chamber, 1-2 way small with solution filter, gamma stable, DEHP-free

SV21ga - MAX21 - C3Sgap - C7S - B1



## GM21gap

Assembled chamber, 1-2 way small, gamma stable, DEHP-free

SV21ga - MAX21 - C3Sgap - B1



## GM21Fmgap

Assembled chamber, 1-2 way small with solution filter, gamma stable, DEHP-free

SV21ga - MAX21 - C3Sgap - MGn - B1 - C7S



## A1Fgap

Assembled chamber, 1 way small with solution filter, gamma stable, DEHP-free

C1Sga - C3Sgap - C7S - B1



## A1gap

Assembled chamber, 1 way small, gamma stable, DEHP-free

C1Sga - C3Sgap - B1



## AB12Fbgap

Assembled chamber, 1-2 way big with solution filter and blue on-off cap, gamma stable, DEHP-free

B2ga - B3b - B4 - B5gap - B6 - B1



## AB12Fmgap

Assembled chamber, 1-2 way big with solution filter and microripper, gamma stable, DEHP-free

B2ga - B3 - B5gap - B6 - B1 - MGn



## AB12gap

Assembled chamber, 1-2 way big, gamma stable, DEHP-free

B2ga - B3 - B4 - B5gap - B6 - B1



## AB12gas

Assembled chamber, 1-2 way big colored, gamma stable, DEHP-free

B2ga - B3 - B5gaps - B6 - B1



## AB12mganp

Assembled chamber, 1-2 way big with microripper, gamma stable, DEHP-free

B2ga - B3 - B4 - B5gap - B6 - B1 - MGn



## AB2Fgap

Assembled chamber, 2 way big with solution filter, gamma stable, DEHP-free

B2ga - B4 - B5gap - B6 - B1



## AB2gap

Assembled chamber, 2 way big, gamma stable, DEHP-free

B2ga - B5gap - B6 - B1

# Assembled chambers



## N11Fgap

Assembled chamber, 1 way big  
with solution filter, gamma stable,  
DEHP-free

B11ga - B4 - B5gap - B1



## N11gap

Assembled chamber, 1 way big,  
gamma stable, DEHP-free

B11ga - B5gap - B1



## AB12Fgap

Assembled chamber, 1-2 way big  
with solution filter, gamma stable,  
DEHP-free

B2ga - B3 - B4 - B5gap - B6 - B1

# Bag Connectors



## PCS2

Silicone stopper 20mm  
Silicone



## SKI

Plug (stopper port), tube I.D.  
6,0mm  
PC



## SKIP

Plug (stopper port), tube I.D.  
6,0mm  
PP



## SKI66

Plug (stopper port), tube I.D.  
6,6mm  
PC



## SKI66P

Plug (stopper port), tube I.D.  
6,6mm  
PP



## SVI

Injection port, tube I.D. 6,0mm  
PC - Rubber - AL



## SPIP

Injection port latex-free, 6,0mm  
PP - Polyisoprene



## SPI66

Injection port latex-free, tube I.D.  
6,6mm  
PC - Polyisoprene



## SPI66P

Injection port latex-free, tube I.D.  
6,6mm  
PP- Polyisoprene



## SPIm

Injection port with membrane,  
tube I.D. 6.0mm  
PC - Polyisoprene



## SV48ct

Injection port with cap, tube I.D.  
4,8mm  
PC - Rubber - AL - PP



## SV59ct

Injection port with cap, EVA or  
PP tube I.D. 6,0mm  
PC - Rubber - AL - PP



# Bag Connectors



## SV66ct

Injection port with cap, tube I.D. 6,6mm

PC - Rubber - AL - PP



## SVIP

Injection port, tube I.D. 6,0mm

PP - Rubber - AL



## SVIct

Injection port with cap, tube I.D. 6,0mm

PC - Rubber - AL -PP



## SVIPct

Injection port with cap, tube I.D. 6,0mm


PP - Rubber - AL



## SVIS2ct

Injection port with cap, tube I.D. 6,0mm


PC - Silicone - AL - PP



## STL15

Long twist-off

PVC



## STL15e

Long twist-off


EVA



## STWS

PVC Twist-off, welded


PVC



## STWSStp3

PVC DEHP-free Twist-off, welded


PVC DEHP-free



## STWSb0

Free PVC twist-off

TPE



## STWS6b0

Free PVC twist-off

TPE



## SPI

Injection port latex-free, tube I.D. 6,0mm

PC - Polyisoprene

# Components



## B5gap

Drip chamber, big, gamma stable, DEHP-free

PVC DOP-free



## C12Sga

2 way spike, small, gamma stable

PVC



## C5S

Air filter, small, blue

PP - Membrane



## C7S

Solution filter small

PE/LD - PA



## MAX21

Air filter with on-off cap, blue

PP - Membrane



## SV21ga

2 way spike, small, gamma stable

PVC



## B1

Protective cap for spike

PE/LD



## B11ga

1 way spike, big, gamma stable

PVC



## B2ga

2 way spike, big, gamma stable

PVC



## B4

Solution filter big

PE/LD - PA



## B6

Air filter, big, red

PP - Membrane




## CC20ga

Cover for drip chamber, gamma stable

PVC

# Components



## MGn

Microdripper adapter, 60 drops/ml

ABS



## IPA

Air filter, small, natural

PP - Membrane



## D4

Protective cap for luer connector D3

PE/LD



## B3

On-off cap, red

PE/LD



## C1Sga

1 way spike, small, gamma stable

PVC



## D3

Luer record connector for tube 3x4,1mm

PVC



## C3Sgap

Drip chamber, small, gamma stable, DEHP-free

PVC DOP-free



# Silicon products



## PCS2

Silicone stopper 20mm  
Silicone



## SVIS2ct

Injection port with cap, tube I.D.  
6,0mm  
PC - Silicone - AL - PP